

# INTENSITY™ Select Combo

INSTRUCTION MANUAL



**This manual is valid for the  
InTENSity™ Select Combo TENS/IF/MIC/EMS Stimulator**

This user manual is published by Richmar.

Richmar does not guarantee its contents and reserves the right to improve and amend it at any time without prior notice. Amendments may however be published in new editions of this manual.

All Right Reserved. Rev. V2 © 2019



United States Federal Law restricts this device to sale by or on the order of a physician or licensed practitioner

**Conformity to safety standards**

Richmar declares that the device complies with following normative documents:

**IEC60601-1, IEC60601-1-2, IEC60601-2- 10, IEC60601-1-4,  
ISO10993-5, ISO10993-10, ISO10993-1**

# TABLE OF CONTENTS

1. GENERAL INFORMATION	4
1.1 General Description	
1.2 Medical Background	
1.3 Indication for Use	
2. SAFETY INFORMATION	8
2.1 Contraindications	
2.2 Warnings, Cautions, Adverse Reactions	
3. PRESENTATION	13
3.1 Front and Rear Panel	
3.2 LCD Display	
4. SPECIFICATION	15
4.1 Accessories	
4.2 Technical Information	
4.3 The Waveforms of the Stimulation Programs	
5. INSTRUCTIONS FOR USE	21
5.1 Battery	
5.2 Connect Electrodes to Lead Wires	
5.3 Connect Lead Wires to Device	
5.4 Electrode	
5.5 Turn ON	
5.6 Select the Therapeutic Program	
5.7 Steps to a New Program	
5.8 Safety Lock Feature	
5.9 Stop the Treatment	
5.10 Turn OFF	
5.11 Low Battery Indicator	
6. PROGRAM	33
7. CLEANING AND CARE	34
7.1 Tips for Skin Care	
7.2 Cleaning the Device	
7.3 Electrodes	
7.4 Cleaning the Electrodes Cords	
7.5 Maintenance	
8. TROUBLESHOOTING	37
9. STORAGE	39
10. DISPOSAL	39
11. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES	40
12. GLOSSARY OF SYMBOLS	44
13. WARRANTY	45

# **1. GENERAL INFORMATION**

## **1.1 General Description**

InTENSity™ Select Combo stimulator is a portable electro-therapy device featuring two therapeutic modes: Transcutaneous Electrical Nerve Stimulation (TENS) and Electrical Muscle Stimulation (EMS), Interferential (IF), and Micro current (MIC) which are used for pain relief and electrical muscle stimulation. The stimulator sends a gentle electrical current to underlying nerves and muscle group via electrodes applied on the skin. The parameters of the device are controlled by the buttons on the front panel. The intensity level is adjustable according to the needs of the individual patient.

## **1.2 Medical Background**

### **EXPLANATION OF PAIN**

Pain is a warning system and the body's method of telling us that something is wrong. Pain is important; without it abnormal conditions may go undetected, causing damage or injury to vital parts of our bodies. Even though pain is a necessary warning signal of trauma or malfunction in the body, nature may have gone too far in its design.

Aside from its value in diagnosis, long-lasting persistent pain serves no useful purpose. Pain does not begin until the coded message travels to the brain where it is decoded, analyzed, and then reacted to. The pain message travels from the injured area along the small nerves leading to the spinal cord. Here the message is switched to different nerves that travel up the spinal cord to the brain. The pain message is then interpreted, referred back and the pain is felt.

### **EXPLANATION OF TENS**

Transcutaneous Electrical Nerve Stimulation (TENS) is a non-invasive, drug free method of controlling pain. TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain. TENS does not work for everyone; however, in most patients it is effective in reducing or eliminating the pain, allowing for a return to normal activity.

## **HOW TENS WORKS**

There is nothing "magic" about Transcutaneous Electrical Nerve Stimulations (TENS). TENS is intended to be used to relieve pain. The TENS unit sends comfortable impulses through the skin that stimulate the nerve (or nerves) in the treatment area. In many cases, this stimulation will greatly reduce or eliminate the pain sensation the patient feels. Pain relief varies by individual patient, mode selected for therapy, and the type of pain. In many patients, the reduction or elimination of pain lasts longer than the actual period of stimulation (sometimes as much as three to four times longer). In others, pain is only modified while stimulation actually occurs. You may wish to discuss this method of pain management treatment with your physician or therapist.

## **EXPLANATION OF EMS**

Electrical Muscle Stimulation (EMS) is an internationally accepted and proven way of treating muscular injuries. It works by sending electronic pulses to the muscle needing treatment; this causes the muscle to exercise passively. It is a product derived from the square waveform, originally invented by John Faraday in 1831. Through the square wave pattern it is able to work directly on muscle motor neurons. This device has low frequency and in conjunction with the square wave pattern allows direct work on muscle groupings. This is widely used in hospitals and sports clinics for the treatment of muscular injuries and for the reeducation of paralyzed muscles, to prevent atrophy in affected muscles and improving muscle tone and blood circulation.

## **HOW EMS WORKS**

The EMS units send comfortable impulses through the skin that stimulate the nerves in the treatment area. When the muscle receives this signal it contracts as if the brain has sent the signal itself. As the signal strength increases, the muscle flexes as in physical exercise. When the pulse ceases, the muscle relaxes and then this cycle is repeated until therapy is completed.

The goal of electrical muscle stimulation is to achieve contractions or vibrations in the muscles. Normal muscular activity is controlled by the central and peripheral nervous systems, which transmit electrical signals to the muscles. EMS works similarly but uses an external source (the stimulator) with electrodes attached to the skin

for transmitting electrical impulses into the body. The impulses stimulate the nerves to send signals to a specifically targeted muscle, which reacts by contracting, just as it does with normal muscular activity.

## **WHAT IS IF?**

Interferential Stimulation IF is an anti-inflammatory based treatment modality. Interferential stimulation is characterized by two alternation-current sine waves or square waves of differing frequencies that “work” together to produce an interferential current that is also known as a beat pulse or alternating modulation frequency. One of the two currents is usually held at 4,000 Hz, and the other can be held constant or varied over a range of 4,001 to 4,100 Hz. Because of the frequency, the interferential wave meets low impedance when crossing the skin to enter deep into soft tissues. The interferential currents reportedly can stimulate sensory, motor, and pain fibers. These large impulse fibers interfere with the transmission of pain messages at the spinal cord level. This deep tissue penetration stimulates parasympathetic nerve fibers for increased blood flow and edema reduction. It utilizes the low electric-current to stimulate muscle nerves to achieve the symptomatic relief of chronic intractable pain, post-traumatic pain, and post-surgical pain.

## **WHAT IS MICROCURRENT?**

Microcurrent stimulation is a type of therapy in which very low current is sent into the cells of the body. Microcurrent is a very faint current that is so small it is measured in millionths of an amp (Microamps). Human cells generate a current that is in the micro amp range which is why you can't feel it - the current is so low it doesn't stimulate the sensory nerves.

Microcurrent is a physiological electric modality that increases ATP (energy) production in the cells of your body. This dramatically increases the tissue's healing rate. The immediate response to the correct microcurrent frequency suggests that other mechanisms are involved as well. The exact effects or changes in the tissue are unmistakable; scars will often suddenly soften; trigger points often become less painful within minutes when the “correct” frequency is applied. In many situations the changes seen seem to be long lasting and in many cases permanent

### 1.3 Indication For Use

InTENSity™ Select Combo TENS/IF/MIC/EMS Stimulator may be used for the following conditions:

#### TENS/IF/MIC:



1. Symptomatic relief of chronic intractable pain.
2. Acute post-traumatic pain.
3. Acute post-surgical pain.

#### EMS:

1. Relaxation of muscle spasm.
2. Increase of blood flow circulation.
3. Prevention of disuse atrophy.
4. Muscle re-education.
5. Maintaining or increasing range of motion.
6. Immediate post-surgical stimulation of lower leg muscles to prevent venous thrombosis.

## 2. SAFETY INFORMATION

Read instruction manual before operating. Be sure to comply with all “Contraindications”, “Warnings”, “Cautions” and “Adverse reactions” in the manual. Failure to follow instructions may cause harm to user or device.

Safety Symbols Used in this Manual	
 WARNING	Indicates a potentially hazardous situation which, if not avoided, could result in serious injury and equipment damage.
 CAUTION	Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the device or other property.

### 2.1 Contraindications

1. This device should not be used for symptomatic local pain relief unless etiology is established or unless a pain syndrome has been diagnosed.
2. This device should not be used when cancerous lesions are present in the treatment area.
3. Stimulation should not be applied over swollen, infected, inflamed areas or skin eruptions (e.g. phlebitis, thrombophlebitis, varicose veins, etc.).
4. Electrodes must not be applied to sites that might cause current/ stimulation to flow through the carotid sinus region (anterior neck) or trans-cerebrally (through the head).
5. **DO NOT** use this device if the patient has a demand-type cardiac pacemaker or any implanted defibrillator.
6. This device should not be used over poorly enervated areas.
7. This device should not be used on patients with epilepsy.
8. This device should not be used on patients with serious arterial circulatory problems in the lower limbs.
9. This device should not be used on patients with abdominal or inguinal hernia.
10. **DO NOT** use this device if you have heart disease without consulting your physician.



## 2.2 Warnings, Cautions and Adverse Reactions



### WARNINGS:

1. This device should be used only under the continued supervision of a licensed physician or practitioner.
2. The long-term effects of chronic electrical stimulation are unknown. Electrical stimulation devices **DO NOT** have any curative value.
3. TENS is a symptomatic treatment and, as such, suppresses the sensation of pain, which would otherwise serve as a protective mechanism.
4. Safety has not been established for the use of therapeutic electrical stimulation during pregnancy. **DO NOT** use during pregnancy unless directed by your physician.
5. Electrical stimulation is not effective for pain of central origin, such as a headache.
6. Electrical monitoring equipment (such as ECG monitors and ECG alarms) may not operate properly when electrical stimulation is in use.
7. Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
8. Stimulation should not be applied over the neck or mouth. Severe spasm of the laryngeal and pharyngeal muscles may occur and the contraction may be strong enough to close the airway or cause difficulty in breathing.
9. Stimulation should not be applied transthoracically. Introduction of electrical current into the heart may cause cardiac arrhythmias.
10. Stimulation should not take place while the user is connected to high-frequency surgical equipment, it may cause burn injuries on the skin under the electrodes, as well as problems with the stimulator.
11. **DO NOT** use the stimulator in the vicinity of shortwave or microwave therapy equipment, since this may affect the output power of the stimulator.
12. **NEVER** use in environments with high humidity such as in the bathroom or when having a bath or shower.
13. Caution should be used in applying electrical stimulation to patients suspected of having heart disease. Further clinical data is needed to show there are no adverse results.
14. **NEVER** use near the heart. Stimulation electrodes should **NEVER** be placed anywhere on the front of the thorax (marked by ribs and breastbone), take extreme caution not to place near or on the two large pectoral muscles. Here it can increase the risk of ventricular fibrillation and lead to cardiac arrest.

15. Electrodes should not be placed over the eyes, in the mouth, near the genitals or internally.
16. **NEVER** use on the areas of the skin which lack normal sensation.
17. Apply the electrodes to clean, dry, and unbroken skin only.
18. Keep electrodes separate during treatment, electrodes in contact with each other could result in improper stimulation or skin burns.
19. Keep the stimulator out of reach of children.
20. Consult your doctor if you have any questions or concerns before using this device.



### **CAUTIONS:**

1. Federal law (USA) restricts this device to sale by or on the order of a physician.
2. This device is for single patient use only.
3. Keep yourself informed of the contraindications.
4. This device is not intended for use on an unattended patient who is non-compliant, emotionally disturbed, has dementia, or low IQ.
5. Read, understand, and practice the warnings, cautions and operating instructions. Know the limitations and hazards associated with using this device. Observe the cautionary and operational decals placed on the unit. Always follow the operating instructions prescribed by your healthcare practitioner.
6. The instruction of use was listed; any improper use may be dangerous.
7. Do not use this device for undiagnosed pain syndromes until consulting a physician.
8. Patients with an implanted electronic device, such as a cardiac pacemaker, implanted defibrillator, or any other metallic or electronic device should not use this device without first consulting a doctor.
9. Stimulation delivered by this device may be sufficient to cause electrocution. Electrical current of this magnitude must not flow through the thorax or across the chest because it may cause cardiac arrhythmia.
10. **DO NOT** place electrodes on the front of the throat as spasm of the Laryngeal and Pharyngeal muscle may occur. Stimulation over the carotid sinus (neck region) may close the airways, make breathing difficult, and may have adverse effects on the heart rhythm or blood pressure.
11. **DO NOT** place electrodes on your head or at any sites that may cause the electrical current to flow transcranially (through the head).
12. Patients with heart disease, epilepsy, cancer or any other health condition should not use this device without first consulting a physician.

13. Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or silicone rubber. If rash develops or pain persists, discontinue use and consult a doctor.
14. Electrode placement and stimulation settings should be based on the guidance of prescribing practitioner.
15. Effectiveness is highly dependent upon the patient and the selection of therapy by a person qualified in the management of pain.
16. Isolated cases of skin irritation have occurred at the site of the electrode placement following long-term application. If this occurs, discontinue use and consult your physician.
17. The electrodes are only to be placed on healthy skin. Avoid skin irritation by ensuring that good contact is achieved between electrodes and skin.
18. If the stimulation levels are uncomfortable or become uncomfortable, reduce the stimulation intensity to a comfortable level and contact your physician if problems persist.
19. This device should not be used while driving, operating machinery, close to water or during any activity in which involuntary muscle contractions may put the user at undue risk for injury.
20. **NEVER** use the device in rooms where aerosols (sprays) are used or pure oxygen is being administered.
21. **DO NOT** use it near any highly flammable substances, gases or explosives.
22. **DO NOT** use this device at the same time as other equipment which sends electrical pulses to your body.
23. **DO NOT** confuse the electrode cables and contacts with your headphones or other devices, and **DO NOT** connect the electrodes to other devices.
24. **DO NOT** use sharp objects such as pencil point or ballpoint pen to operate the buttons on the control panel.
25. Inspect applicator cables and associated connectors before each use.
26. Turn the device off before applying or removing the electrodes.
27. Electrical stimulator should be used only with the leads and electrodes recommended for use by the manufacturer.
28. This device has no AP/APG protection. **DO NOT** use it in the presence of explosive atmosphere of flammable mixture.

## **ADVERSE REACTIONS:**

1. Skin irritation from the electrode gel and electrode burns are potential adverse reactions. If skin irritation occurs, discontinue use and consult your physician.

**Note:** Always use electrodes that are legally marketed and sold in the United States under 510K guidelines.

2. If the stimulation levels are uncomfortable, reduce the stimulation intensity to a comfortable level and contact your physician if any problems persist.


## 3. PRESENTATION

### 3.1 Front and Rear Panel



1. Output socket: electric signal output after connection of the cable with adhesive electrodes on channel 1.
2. Output socket: electric signal output after connection of the cable with adhesive electrodes on channel 2.
3. Increases [ $\blacktriangle$ ] the output intensity of channel 1. It is also used to set the application program and the parameter of the waveform in the setting state.
4. Decreases [ $\blacktriangledown$ ] the output intensity of channel 1. It is also used to set the application program and the parameter of the waveform in the setting state, and to unlock the current treatment program.
5. Therapeutic mode [**M**] selection button; it also stops treatment and exits you out of the setting mode to return to the user interface.
6. LCD display: Shows the operating state of the device.
7. Increases [ $\blacktriangle$ ] the output intensity of channel 2. It is also used to set the application program and the parameter of the waveform in the setting state.
8. Decreases [ $\blacktriangledown$ ] the output intensity of channel 2. It is also used to set the application program and the parameter of the waveform in the setting state, and unlocks the current treatment program.
9. Parameter Selection [**S**] button: press the button to enter the setting state; you can select the different parameters in conjunction with the up [ $\blacktriangle$ ] and down [ $\blacktriangledown$ ] buttons.
10. Press [ $\text{⏻}$ ] to turn ON the device. Press [ $\text{⏻}$ ] and hold for approximately 3 seconds to turn off the device.
11. Belt Clip.
12. Battery compartment cover.
13. AC Adapter connecting port

### 3.2 LCD Display



1. Displays therapeutic mode.
2. Displays therapeutic program for TENS and EMS therapeutic mode
3. Displays therapeutic program for IF and MIC therapeutic mode or displays the cycle time for TENS, IF and MIC therapeutic mode in setting state.
4. Timer symbol
5. EMS waveform of ramp up and ramp down time.
6. Display of waveform pulse width.
7. Displays channel 1.
8. Displays the output intensity for channel 1 (CH1) Displays waveform pulse width or EMS waveform of contraction (working) time in setting state.
9. EMS waveform of contraction (working) time.
10. The device is locked indicator
11. Displays the waveform pulse rate.
12. Displays the output intensity for channel 2 (CH2). Displays waveform pulse rate or EMS waveform of relaxation time in setting state.
13. Displays channel 2.
14. Displays the EMS waveform of relaxation time.
15. Displays the treatment time or EMSE waveform of ramp up and ramp down time.
16. Low-battery indicator.

## 4. SPECIFICATIONS

### 4.1 Accessories

NO	DESCRIPTION	QTY
1	TENS/IF/MIC/EMS Device (Item: DI8195)	1 each
2	Pair of Lead wires (Item: WW3005)	2/pk
3	2" x 2" White Cloth electrodes (Item: EP2020WC2-INTM)	4/pk
4	9V TENScell alkaline battery (Item: TA5013-I)	1 each
5	Instruction manual	1 each
6	Carrying case (Item: CC5082)	1 each
7	Wall (AC) adapter (Item: DI0001X)	1 each

### 4.2 Technical Information

Channel	Dual, isolated between channels
Power supply	9.0 V Alkaline battery, type: 6LR61 Adapter output:9.0Vdc, 800mA (optional)
Operating Conditions	5°C to 40°C (41°F to 104°F) with a relative humidity of 30%-75%,atmospheric pressure from 700 to 1060 Hpa
Storage Conditions	-10°C to 50°C (14°F to 122°F) with a relative humidity of 10%-90%,atmospheric pressure from 700 to 1060 Hpa
Dimensions	4.5 × 2.55 × 0.9 inches (L*W*H)
Weight	0.28 lbs.(With battery)
Tolerance	There may be a ±5% tolerance of all setting and ±10% tolerance of output of intensity.
Timer	Adjustable from 1 to 60 minutes or continuous. Adjusts in 1 minute steps. Treatment time countdown is automatic.
Electrode Detection Function	The amplitude level will be reset to 0mA when the amplitude level is 12mA or greater and if an open circuit at either channel is detected.

## Technical Specifications for Transcutaneous Electrical Nerve Stimulator (TENS) Mode

Waveform	Mono-phase square pulse wave
Pulse Amplitude	Adjustable, 0 - 105mA at 1000 ohm Load each channel, 1mA/Step.
Pulse Width	From 50 to 300 $\mu$ s, 10 $\mu$ s/step
Pulse Rate	Adjustable from 1 to 150 Hz, 1 Hz/step
Burst (B)	Burst rate: Adjustable 0.5-5Hz, Pulse Width Adjustable 50-300 $\mu$ s Frequency Fixed = 100 Hz
Normal (N)	The pulse rate and pulse width are adjustable. It generates continuous stimulation based on the setting value.
Pulse Width Modulation (M)	The pulse width is automatically varied in a cycle time. The pulse width is decreased from its original setting to 60% in setting cycle time, and then increased from 60% to its original setting in next setting cycle time. In this program, pulse rate (1 to 150Hz), pulse width (50 to 300 $\mu$ s) and cycle time (5 to 30 seconds) are fully adjustable.
Pulse Width Modulation (M1)	The pulse rate is automatically varied in a cycle time. The pulse rate is decreased from its original setting to 60% in setting cycle time, and then increased from 60% to its original setting in next setting cycle time. In this program, pulse rate (1 to 150Hz), pulse width (50 to 300 $\mu$ s) and cycle time (5 to 30 seconds) are fully adjustable.



## Technical Specifications for Electrical Muscle Stimulation (EMS) Mode

Waveform	Mono-phase square pulse
Pulse Amplitude	Adjustable 0-105mA at 1000 ohm load each channel, 1 mA/step.
Pulse Width	Adjustable from 50 to 300 $\mu$ s, 10 $\mu$ s/step
Pulse Rate	Adjustable from 1 to 150 Hz, 1 Hz/step
Contraction time	Adjustable 1-60 seconds, 1 sec./step
Relaxation (OFF)	Adjustable 0-60 seconds, 1 sec./step
Ramp time	Adjustable, 1-6 seconds, 1 sec./step, The "On" time will increase and decrease in the setting value.
Synchronous (s)	Stimulation of both channels occurs synchronously. The "ON" time including "Contraction", "Ramp Up" and "Ramp Down" time. ON TIME=Contraction + Ramp up + Ramp down.
Alternate (A)	The Stimulation of CH2 will occur after the 1st operation of CH1 is completed. In this program, The "ON" time includes "Contraction", "Ramp Up" and "Ramp Down" time. The OFF Time should be equal to or more than the ON time. ON TIME=Contraction + Ramp up + Ramp down OFF TIME $\geq$ ON TIME
Delay (D)	The Stimulation of CH2 will occur after the 1st operation of CH1 is started + Delay Time. In this program, the "ON" time includes "Contraction", "Ramp Up" and "Ramp Down" time. The OFF Time should be equal to or more than the ON Time + Delay Time. ON TIME=-Contraction + Ramp up + Ramp down.

## Technical Specifications for Interferential (IF) mode


Waveform	Bi-phase square pulse
Pulse Amplitude	Adjustable 0-70mA peak to peak at 1000 ohm load
Pulse Width	125µs
Pulse Rate	Channel 1 -Fundamental frequency: 4000 Hz fixed Channel 2 -Selectable frequency: 4001 to 4150 Hz Interference frequency:1 to 150 Hz

## Technical Specifications for Microcurrent (MICRO) mode


Waveform	Mono-phase square pulse
Pulse Amplitude	Adjustable 0-70mA peak at 1000 ohm load each channel, 0.01mA/Step
Pulse Width (P.W.)	Adjustable, from 2 to 200µs, 1µs/step $P.W. \leq \frac{1}{2 \times P.R.}$
Pulse Rate (P.R.)	Adjustable, from 1-150 Hz/step $P.W. \leq \frac{1}{2 \times P.R.}$
Constant (P1)	Constant stimulation based on setting value. Only pulse width, pulse rate and timer are adjustable in this program. "Constant" is equal to the "Normal" mode of a TENS therapeutic mode
Pulse Width Modulation (P2)	The pulse width is automatically varied in a cycle time. The pulse width is decreased from its original setting to 60% in setting cycle time, and then increased from 60% to its original setting in nest setting cycle time. In this program, pulse rate (1 to 150Hz), pulse width (2 to 200µs) and cycle time (5 to 30 sec) are fully adjustable.
Pulse Rate Modulation (P3)	The pulse rate is automatically varied in a cycle time. The pulse rate is decreased from its original setting to 60% in setting cycle time, and then increased from 60% to its original setting in nest setting cycle time. In this program, pulse rate (1 to 150Hz), pulse width (2 to 200ms) and cycle time (5 to 30 sec) are fully adjustable.

## 4.3 The Waveforms of the Stimulation Programs


### Burst (B)




### Normal (N)




### Pulse Width Modulation



### Pulse Rate Modulation




### Synchronous (S)




### Alternate (A)



### Delay (D)



### Interferential



### Microcurrent (Constant)




## 5. INSTRUCTIONS FOR USE

### 5.1 Battery

#### 5.1.1 Check/Replace the Battery

Over time, in order to ensure the functional safety of the device, changing the battery is necessary.

1. Slide the battery compartment cover to open.
2. Insert the 9V battery into the battery compartment.
3. Make sure you are installing the battery properly. Be sure to match the positive and negative ends of the battery to the marking in the battery compartment of the device.
4. Press and push down the battery following the direction of the arrow indicated on the photo.
5. Replace the battery compartment cover and press to close.
6. If battery requires replacement, slide the battery compartment cover to open. Pull up the battery following the directions of the arrow indicated on the photo and insert a new 9V battery according to the above steps 2 to 5.



#### 5.1.2 Disposal of Battery

Dispose of used batteries according to the current federal, state and local regulations. As a consumer, you are obligated by law to discard spent batteries appropriately.





### **CAUTION FOR BATTERIES:**

1. Swallowing a battery may be fatal. Keep the battery and the device out of the reach of children. If a battery is swallowed, consult a physician immediately.
2. If a battery has leaked, avoid contact with skin, eyes and mucus membranes. Rinse the affected areas with clear water immediately and contact a physician immediately.
3. Battery should not be charged, dismantled, thrown into fire or short-circuited.
4. Protect battery from excess heat.
5. Remove batteries from the unit if they are depleted or if you are not using the unit for prolonged periods of time. This prevents damage caused by leaking battery.
6. Always replace with the same type battery.

### **5.2 Connect Electrode to Lead Wires**

Insert the lead wire connector into electrode connector (standard 0.08 inch female connection). Make sure the connectors are completely pushed together showing no exposed metal of the pins.



### **CAUTION:**

Always use the electrodes with CE mark, or which are legally marketed in the U.S. under 510(K) procedure.

## 5.3 Connect Lead Wires to Device

1. Before proceeding to this step, be sure the device is completely turned OFF.
2. Insert the wires provided with the system into the jack sockets located on top of the device.
3. Holding the insulated portion of the connector, push the plug end of the wire into one of the jacks (see drawing); one or two sets of wires may be used.
4. This device has two output receptacles controlled by Channel 1 and Channel 2 at the top of the unit. You may choose to use one channel with one pair of lead wires or both channels with two pairs of lead wires. Using both channels gives the user the advantage of stimulating two different areas at the same time.



### CAUTION:

**DO NOT** insert the plug of the patient lead wire into any AC power supply socket.


## 5.4 Electrode

### 5.4.1 Electrode Options

The electrodes are disposable and should be routinely replaced before they start to lose their adhesive nature. Replacement electrodes should be re-ordered through or on the advice of your physician to ensure proper quality. Follow application procedures outlined in electrode package to maintain optimal stimulation and to prevent skin irritation.

### 5.4.2 Place Electrode on Skin

Apply electrodes to the exact site indicated by your physician or therapist. Before applying electrodes, be sure the skin surface over which electrodes are placed is thoroughly cleaned and dry.



Make sure the electrodes are pressed firmly to the skin and make good contact between the skin and the electrodes. Place the electrodes over the skin; attach them properly, firmly and evenly. Electrodes should be placed at least 2" but no more than 6" apart, per channel.

#### CAUTION FOR ELECTRODES:

1. Before applying the self-adhesive electrodes, it is recommended that you wash, degrease and dry the skin first.
2. **DO NOT** turn on the device when the self-adhesive electrodes are not positioned on the body.
3. **NEVER** remove the self-adhesive electrodes from the skin while the device is turned on. You will feel an uncomfortable electrical shock.
4. It is recommended that, at a minimum, 1.5" x 1.5" self-adhering, square electrodes are used at the treatment area.

### 5.4.3 Electrode Placement


The placement of electrodes can be one of the most important parameters in achieving success with this therapy. Of utmost importance is the willingness of the physician to try the various styles of electrode placement to find which method best fits the needs of the individual patient.

Every patient responds to electrical stimulation differently and their needs may vary from the conventional settings suggested here. If the initial results are not positive, speak to your physician about alternative stimulation settings and/or electrode placements. Once an acceptable location has been achieved, mark down the electrode sites and the device settings, so the patient can easily continue treatment on their own.



## 5.5 Turning On the Device

Before using the device for the first time, you are strongly advised to take careful note of the contraindications and safety measures detailed at the beginning of this manual (pages 7-11 and throughout manual), as this powerful equipment is neither a toy nor a gadget!

To turn on the device, PRESS and RELEASE the ON/OFF  button. The operation page will appear on the LCD screen.

## 5.6 Select the Therapeutic Part Program


There are four therapeutic modes available -TENS, EMS, IF, and MIC. The therapeutic mode can be selected by pressing the [M] button control.

 **CAUTION:**



Consult your physician for your suitable therapeutic mode.

## 5.7 Steps to Set a New Program

### 5.7.1 TENS Setting


<p><b>Step 1: Set the Therapeutic Program</b></p> 	<p>There are 4 programs available in TENS therapeutic mode available - Burst (B), Normal (N), Pulse Width Modulation (M) and Pulse Rate Modulation (M1). Press [<b>M</b>] until there is a box around TENS, press [<b>S</b>] and “B” will be outlined by a flashing box. To change to either N, M, or M1, press one of the [<b>▲</b>] until there is a flashing box around the desired program.</p>
<p><b>Step 2: Set Timer</b></p>	<p>Press the [<b>S</b>] button to cycle and to enter this setting. The treatment time is adjustable from 1 to 60 minutes or Continuous. Press the up [<b>▲</b>] or down [<b>▼</b>] button controls to adjust the setting. You can set the timer to “Continuous” mode by pressing the up [<b>▲</b>] control when it shows 60 minutes. The output will shut off when time is up. When in Continuous mode, you will need to manually turn off the device.</p>
<p><b>Step 3: Set Pulse Width</b></p>	<p>Pulse Width is adjustable from 50<math>\mu</math>s to 300<math>\mu</math>s. Press the [<b>S</b>] button to enter this menu, then press the up [<b>▲</b>] or down [<b>▼</b>] button to adjust the setting.</p>
<p><b>Step 4: Set Pulse Rate</b></p>	<p>Pulse rate is adjustable from 1 Hz to 150 Hz (0.5 Hz to 5 Hz for Burst). Press the [<b>S</b>] button to cycle and to enter this menu, and then press the up [<b>▲</b>] or down [<b>▼</b>] button to adjust the setting.</p>
<p><b>Step 5: Start Treatment</b></p>	<p>Once you have pressed the [<b>S</b>] button to bring you back to the main screen, you will see a 0mA in the lower right and left corner of the screen. Use the corresponding channels up [<b>▲</b>] or down [<b>▼</b>] to increase the intensity to a comfortable stimulation. If the intensity is uncomfortable at any time press up [<b>▲</b>] or down [<b>▼</b>].</p>


## 5.7.2 EMS Setting

<p><b>Step 1: Set the Therapeutic Program</b></p> 	<p>There are 3 programs in EMS therapeutic mode available: Synchronous (S), Alternate (A) and Delayed (D). Press [M] until there is a box around EMS, press [S] and “S” will be outlined by a flashing box. To change to either A or D, press one of the [▲] until there is a flashing box around the desired program.</p>
<p><b>Step 2: Set Timer</b></p>	<p>Press the [S] button to cycle and to enter this setting. To adjust the timer, press the up [▲] or down [▼] button to choose between 1-60 mins. To choose “Continuous” press [▲] once after 60 is displayed. When in Continuous Timer, you will need to manually turn the device off to stop treatment. Otherwise, device will shut off when time is up.</p>
<p><b>Step 3: Set Pulse Width</b></p>	<p>Pulse Width is adjustable from 50<math>\mu</math>s to 300<math>\mu</math>s. Press the [S] button to enter this menu, then press the up [▲] or down [▼] button to adjust the setting.</p>
<p><b>Step 4: Set Pulse Rate</b></p>	<p>Pulse rate is adjustable from 1 Hz to 150 Hz. Press the [S] button to cycle and to enter this menu, and then press the up [▲] or down [▼] button to adjust the setting.</p>
<p><b>Step 5: Set Delay Time (Optional)</b></p> 	<p>Only the “Delay” therapeutic program has this parameter setting. Press the [S] button to cycle and to enter this menu, and then press the up [▲] and down [▼] button to adjust the setting.</p>
<p><b>Step 6: Set Ramp Time</b></p>	<p>The ramp time controls the time of the output current that increases from 0 to the setting level; and from the setting value back to 0. The ramp time is adjustable from 1-6 sec. Press the [S] button until “Ramp” flashes on the screen. Then press up [▲] and down [▼] button to desired time.</p>

<p><b>Step 7: Set Contract Time</b></p>	<p>The contract time controls the time of stimulation and can be adjusted by pressing the [S] button cycle to enter this menu, and then pressing the up [▲] and down [▼] button to adjust the setting. The stimulation of Channel 1 and Channel 2 is cycled on and off by the contraction and relaxation settings. The range is adjustable from 1 to 60 seconds.</p> <p><b>NOTE:</b> Contract time does not include the ramp up and ramp down time. ON time = Ramp up + Contract time + Ramp down.</p>
<p><b>Step 8: Set Relaxation (OFF) time</b></p>	<p>The Off Time controls the time of relaxation and can be adjusted. Press the [S] button to cycle and to enter this menu, and then press the up [▲] and down [▼] button to adjust the setting. The stimulation of Channel 1 and Channel 2 is cycled on and off by the contraction and relaxation settings. The range is adjustable from 1 to 60 seconds.</p> <p><b>NOTE:</b> In Alternate (A) program, the OFF Time should be equal to or more than the ON Time. (OFF TIME ≥ ON TIME)</p>
<p><b>Step 9: Start Treatment</b></p>	<p>Once you have pressed the [S] button to bring you back to the main screen, you will see a 0mA in the lower right and left corner of the screen. Use the corresponding channels up [▲] or down [▼] to increase the intensity to a comfortable stimulation. If the intensity is uncomfortable at any time press down [▼] to decrease intensity.</p>

### 5.7.3 IF Setting

<p><b>Step 1: Set the Therapeutic Program</b></p>  <p>The screenshot shows a digital display with the following information: At the top, 'TENS' is selected with a box around it, followed by 'IF', 'MIC', and 'EMS'. Below this, 'P4' is displayed with a small circular icon to its left. In the center, '80 min.' is shown in a large font. Below that, '0 mA' is displayed on the left and '80' with a small 'tr' superscript on the right.</p>	<p>There are 4 programs in IF therapeutic mode available. Press [M] until there is a box around IF. To choose the specific program under IF, press [S] and "P1" will be displayed on the screen. To change to P2-P4, press one of the [▲] until that program is displayed on the screen.</p>
---	--

<p><b>Step 2: Set Timer</b></p>	<p>Press the [<b>S</b>] button to cycle and to enter this setting. To adjust the timer, press the [<b>▲</b>] and [<b>▼</b>] buttons to choose between 1-60 mins. To choose “Continuous” press [<b>▲</b>] once after 60 is displayed. When in Continuous Timer, you will need to manually turn the device off to stop treatment. Otherwise, device will shut off when time is up.</p>
<p><b>Step 3: Set Cycle Timer</b></p>	<p>Program P4 is the only place you can adjust this setting. Press [<b>S</b>] until the Hz parameter in the lower right corner of the screen is flashing. Then press [<b>▲</b>] and [<b>▼</b>] to adjust to the desired setting.</p>
<p><b>Step 4: Set Interference frequency (optional)</b></p>	<p>Channel 1 has 4000 Hz fixed Fundamental frequency. Channel has selectable frequency from 4001 to 4150Hz; Interference frequency is adjustable from 1 Hz to 150 Hz. Only “P4” has this parameter setting. Press [<b>S</b>] button cycle to enter this menu, and then press the [<b>▲</b>] and [<b>▼</b>] buttons to adjusting the setting.</p>
<p><b>Step 5: Start Treatment</b></p>	<p>Once you have pressed the [<b>S</b>] button to bring you back to the main screen, you will see a 0mA in the lower right and left corner of the screen. Use the corresponding channels up [<b>▲</b>] or down [<b>▼</b>] to increase the intensity to a comfortable stimulation. If the intensity is uncomfortable at any time press down [<b>▼</b>] to decrease intensity.</p>
	<p><b>Note:</b> For IF setting, if using both channels, the electrodes must be criss-crossed in an X pattern as shown to the left for treatment.</p>

## 5.7.4 Microcurrent Setting

<b>Step 1: Set the Therapeutic Program</b>	
	<p>There are 3 Programs in MIC mode: P1-P3 Press [M] until there is a box around MIC. To choose the the specific program under MIC, press [S] and “P1” will be displayed on the screen. To change to P2-P3, press one of the [▲] until that program is displayed on the screen.</p>
<b>Step 2: Set Timer</b>	<p>Press the [S] button to cycle and to enter this setting. To adjust the timer, press the [▲] and [▼] to choose between 1-60 mins. To choose “Continuous” press [▲] once after 60 is displayed. When in Continuous Timer, you will need to manually turn the device off to stop treatment. Otherwise, device will shut off when time is up.</p>
<b>Step 3: Set Pulse Width</b>	<p>Pulse Width is adjustable from 2<math>\mu</math>s to 200<math>\mu</math>s. Press [S] button cycle to enter this menu, then press [▲] and [▼] button to adjust the setting.</p>
<b>Step 4: Set Pulse Rate</b>	<p>Pulse rate is adjustable from 1 Hz to 150 Hz. Press the [S] button to cycle and to enter this menu, and then press the up [▲] and down [▼] button to adjust the setting.</p>
<b>Step 5: Set Cycle Time (Optional)</b>	
	<p>Cycle Time is only adjustable in P2-P3. Press [S] button cycle to enter this menu, and then press the [▲] and [▼] button to adjusting the setting.</p>
<b>Step 6: Start Treatment</b>	<p>Once the [S] button is pressed and all parameters are set, the device will go back to the home screen and display 0.00 mA in the lower left and right corner of the screen. Use the [▲] and [▼] of the corresponding channel to increase the intensity to a comfortable stimulation or the setting recommended by your physician. <b>Note:</b> Most people do not feel Microcurrent Stimulation, as it is administered using Microamps, instead of Milliamps.</p>

 **CAUTION:**

1. If the stimulation levels are uncomfortable or become uncomfortable, reduce the stimulation intensity to a comfortable level and contact your medical practitioner if problems persist.
2. When using TENS, EMS, or IF therapeutic mode, if the electrodes are not placed firmly on the skin or the device is not connected to the electrodes and the stimulator's output intensity surpasses 12mA, the intensity will automatically reset to 0mA.

## 5.8 Safety Lock Feature

The Safety Lock Feature automatically activates by locking out your ability to operate the buttons when you have stopped setting the program or adjusting the intensity for 30 seconds. This is a safety feature to prevent accidental changes to your settings and to prevent accidental increases to the intensity levels. You can press either one of the down [▼] buttons to unlock the device.

## 5.9 Stop the Treatment

When you have activated the treatment timer, you can press the [M] button or the down [▼] button to control or stop the treatment.

 **CAUTION:**

Default state, if the button is locked, you can press only one of the down [▼] buttons to unlock, and then press the [M] button or the down [▼] button to control stop the treatment.

## 5.10 Turn Off

Press and HOLD the [⏻] button for approximately 3 seconds to turn OFF the device. **Note:** If device is locked [m-O] and treatment is active, in order to turn off the device, you must press the [▼] button before pressing the [⏻] to turn device off.

 **CAUTION:**

1. If there is no operation in the panel for 2 minutes in the waiting state, the device will turn off automatically.
2. In shutdown state, keep pressing the channel 2 down [▼] button first, and then press the [⏻] button at the same time to restore factory parameter settings.

### **5.11 Low Battery Indicator**

A battery symbol is shown flashing on the display when the battery is almost empty. The battery should be replaced with a new battery as soon as possible. However, the unit may continue to operate for an extended period of time depending on the setting and intensity level.



## 6. PROGRAM

Mode	Program	Modulation Method	Frequency	Pulse Width	Treatment Time
TENS	B	Burst	0.5-5Hz	50-300µs	1-60 min, continuous
	N	Continuous	1-150Hz	50-300µs	1-60 min, continuous
	M	Pulse width	1-150Hz	50-300µs	1-60 min, continuous
	M1	Frequency Modulation	1-150Hz	50-300µs	1-60 min, continuous
EMS	S	Synchronous mode	1-150Hz	50-300µs	1-60 min, continuous
	A	Asynchronous mode	1-150Hz	50-300µs	1-60 min, continuous
	D	Delay Mode	1-150Hz	50-300µs	1-60 min, continuous
IF	P1	Frequency Modulation	4kHz 4001-4010Hz	125µs	1-60 min, continuous
	P2	Frequency Modulation	4kHz 4001-4010Hz	125µs	1-60 min, continuous
	P3	Frequency Modulation	4kHz 4080-4010Hz	125µs	1-60 min, continuous
	P4	Frequency Modulation	4kHz 4001-4010Hz	125µs	1-60 min, continuous
MIC	P1	Continuous	1-150Hz	2-200µs	1-60 min, continuous
	P2	Pulse width modulation	1-150Hz	2-200µs	1-60 min, continuous
	P3	Frequency Modulation	1-150Hz	2-200µs	1-60 min, continuous

## 7. CLEANING AND CARE

### 7.1 Tips for Skin Care

Follow these suggestions to avoid skin irritation, especially if you have sensitive skin:

1. Wash the area of skin you will be placing the electrodes on with soap. Rinse thoroughly and dry the area completely before and after placing electrodes.
2. Excess hair may be clipped with scissors; **DO NOT** shave stimulation area.
3. Wipe the area with the skin preparation your clinician has recommended. Let this dry. Apply electrodes as directed.
4. Many skin problems arise from the “pulling stress” from adhesive patches that are excessively stretched across the skin during application. To prevent this, apply electrodes from center outward; avoid stretching over the skin.
5. To minimize “pulling stress”, tape extra lengths of lead wires to the skin in a loop to prevent tugging on electrodes.
6. When removing electrodes, always remove by pulling in the direction of hair growth.
7. It may be helpful to rub skin lotion on electrode placement area during treatment down time when you are not wearing electrodes.
8. **NEVER** apply electrodes over irritated or broken skin.

### 7.2 Cleaning the Device

1. Remove the battery from the device before you clean the device.
2. Clean the device after use with a soft, slightly moistened cloth. For hard to clean situations, you can also moisten the cloth with mild soapy water.
3. **DO NOT** use any chemical cleaners or abrasive agents for cleaning.

## 7.3 Electrodes

1. Use this device only with the leads and electrodes provided by the manufacturer. Use only the electrode placements and stimulation settings prescribed by your physician or therapist.
2. It is recommended, at minimum, 1.5" x 1.5" self-adhering electrodes be used at the treatment area.
3. Inspect your electrodes before every use. Replace electrodes as needed. Reusable electrodes can cause slight skin irritation, lose adhesion properties and deliver less stimulation if overused.

### **TO REMOVE YOUR ELECTRODES:**

1. Lift the corner of the electrode and gently remove it from the skin.
2. It may be helpful to improve repeated electrode application by spreading a few drops of cold water over the adhesive side and turn the surface up to air dry. Over saturation with water will reduce the adhesive properties. This can only be done once and then electrodes need to be replaced.
3. Between uses, place the electrodes back onto the protective sleeve and insert them into the resealable bag and store in a cool dry place.



## **CAUTION:**

1. **DO NOT** pull on the electrode wire. Doing so may damage the wire and electrode.
2. **DO NOT** apply to broken skin.
3. The electrodes should be discarded when they are no longer adhering to the skin.
4. The electrodes are intended for single patient use only.
5. If irritation occurs, discontinue use and consult your clinician.
6. Read the instructions for use of self-adhesive electrodes before application.
7. Always use the electrodes with the requirements of the IEC/EN60601-1, ISO10993-1/-5/-10 and IEC/EN60601-1-2, such as with CE mark, or are legally marketed in the U.S. under 510(K) procedures.

### **7.4 Cleaning the Electrode's Cords**

Clean the electrode cords by wiping them with a damp cloth. Coating them lightly with talcum powder will reduce tangles and prolong their life.

### **7.5 Maintenance**

1. Maintenance and all repairs should only be carried out by an authorized agency. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized persons.
2. The user must not attempt any repairs to the device or accessories. Please contact the retailer for repair.
3. Opening of the equipment by unauthorized agencies is not allowed and will terminate any claim to warranty.
4. Check the unit before each use for signs of wear and/or damage. Replace worn items as required.

## 8. TROUBLESHOOTING

If your device does not seem to be operating correctly, refer to the chart below to determine what may be wrong. Should none of these measures correct the problem, the device should be serviced.

Problem	Possible Cause	Solution
Display fails to light up	Battery contact failure.	Try fresh batteries.
		Ensure batteries are inserted correctly.
		Check contacts are in place.
		Check contacts are not broken.
Stimulation weak or cannot feel any stimulation	Electrodes are dried out or contaminated.	Replace and reconnect.
	Electrode placement.	
	Lead wires old/worn/damaged.	Replace
Stimulation is uncomfortable	Intensity is too high.	Decrease intensity.
	Electrodes are too close together.	Reposition the electrodes.
	Damaged or worn electrodes or lead wires.	Replace.
	Electrode active area size is too small.	Replace electrode with ones that have an active area no less than 16.0cm <sup>2</sup> (4cm x 4cm).
	May not operate the device according to the manual.	Please check the manual before use.
Intermittent output	Lead wires.	Verify connection is secure and firmly seated and no metal pins are exposed.
		Turn down the intensity. Rotate lead wires in socket 90°. If still intermittent, replace lead wire.
		If still intermittent after replacing the lead wire, a component may have failed. Call the your distributor or the manufacturer
	Program option in use	Some programs will seem intermittent. This is expected. Refer to the program on page 33.

Stimulation is ineffective	Improper electrode and applicator placement unknown.	Reposition electrode and applicator.
		Contact clinician.
The skin becomes red and/or you feel a stabbing pain	Using the electrodes on the same site every time.	Reposition the electrodes. If at any time you feel pain or discomfort stop use immediately.
	The electrodes are not sticking onto the skin properly.	Ensure the electrode is stuck securely on the skin.
	The electrodes are dirty.	Replace with new electrodes
	The surface of the electrode was scratched.	Replace with new electrodes.
Output current stops during therapy	The electrode pads come off the skin.	Turn off the device and stick the electrode pad firmly to the skin. If that does not work, replace with new electrodes.
	The cable is disconnected.	Turn off the device and connect the cable.
	The power of the batteries has been exhausted.	Replace them with new batteries.

## 9. STORAGE

1. For prolonged pauses in treatment, store the device in a cool dry room and protect it against heat, sunshine and moisture and remove the battery to avoid battery leaking.
2. Store the device in a cool, well-ventilated place.
3. **NEVER** place any heavy objects on the device.

## 10. DISPOSAL

Used fully discharged batteries must be disposed of in a specially labeled collection container, at a toxic waste collection point or through an electrical retailer. Please dispose of the device in accordance with the laws in your area.



## 11. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and manufacturer's declaration - electromagnetic emissions		
<p><b>The device is intended for use in the electromagnetic environment specified below.</b>  <b>The customer or the user assures that it is used in such an environment.</b></p>		
Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group B	The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	




**Guidance and manufacturer's declaration - electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user should assure that it is used in such an environment.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact	±6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
	±8 kV air	±8 kV air	
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply	Not applicable	Main power quality should be that of a commercial or hospital environment.
	±1 kV for input/output lines		
Surge IEC 61000-4-5	±1 kV differential mode		
	±2 kV common mode		
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle		
	40% UT (60% dip in UT) for 5 Cycles		
	70% UT (30% dip in UT) for 25 Cycles		
	<5% UT (>95% dip in UT) for 5 sec		
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Immunity test	IEC 60501 test level	Compliance level	Electromagnetic environment guidance
			Portable and mobile RF Communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = \left[ \frac{3.5}{\sqrt{1}} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[ \frac{3.5}{E1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[ \frac{7}{E1} \right] \sqrt{P}$ 800 MHz to 2.5 MHz
			<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the Transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range b.</p> 

NOTE 1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

1. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

2. Over the frequency range 150 kHz to 80 MHz, field strength should be less than (Vi) W/m.

**Recommended separation distances between portable and mobile RF communications equipment and the device.**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) as recommended below, according to the maximum output power of the communications equipment








Rate maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \left[ \frac{3.5}{\sqrt{1}} \right] \sqrt{P}$	$d = \left[ \frac{3.5}{\sqrt{1}} \right] \sqrt{P}$	$d = \left[ \frac{3.5}{\sqrt{1}} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## 12. GLOSSARY OF SYMBOLS

	Batch code
	Serial number
	Attention: Read the operating instruction before use! Equipment capable of delivering output values in excess of 10mA r.m.s. or 10V r.m.s. averaged over any period of 5s.
	Electrical devices are recyclable material and should not be disposed of with household waste after their useful life! Help us to protect the environment and save resources and take this device to the appropriate collection points. Please contact the organization which is responsible for waste disposal in your area if you have any questions.
	Type BF Applied Part
	Type of protection against electric shock: Class II Equipment
	Refer to instruction manual

## 13. WARRANTY

**Please contact your dealer in case of a claim under the warranty. If you have to send the unit back to your provider, enclose a copy of your receipt and state what the defect is.**

**The following warranty terms apply:**

1. The warranty period for device is one year from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Repairs or replacement under warranty **DO NOT** extend the warranty period either for the device or for the replacement parts.
3. The following is excluded under the warranty:
  - All damage which has arisen due to improper treatment, e.g. non-observance of the user instruction.
  - All damage which is due to repairs or tampering by the customer or unauthorized third parties.
  - Damage which has arisen during transport from the manufacturer to the consumer or during transport to the retailer.
  - Accessories which are subject to normal wear and tear.

Liability for direct or indirect consequential losses caused by the unit is excluded even if the damage to the unit is accepted as a warranty claim.

1. All products must be returned in original packaging and must contain all components, accessories and user manuals. If any components are missing, you will be responsible for the cost of the replacement component and the 25% restocking fee.
2. All returns must be approved with a Return Authorization Number. Please call our Customer Service Team at (800) 376-7263 to obtain a Return Authorization Number. Provide the following information when calling:
  - Item Number
  - Original Order Number
  - Product Serial Number/Lot Number
  - Reason for Return
3. The Return Authorization Number must be marked clearly on the returned carton and is valid for 10 business days from the date of issue.
4. Returned merchandise must be in the same unit of measure as originally purchased.

5. Return Labels or Call Tags can be issued by our customer service department to return merchandise.
6. Associated fees and return freight charges will apply. All returns of dropshipped items are subject to a restocking fee as well as inbound and outbound freight charges.
7. Returns will not be accepted on items that are:
  - Missing their serial number
  - Special order items
  - Returned more than 30 days after delivery
  - Returned without notification







Manufactured for:



**Richmar**<sup>®</sup>

Compass Health Brands Corp.  
Toll Free 1,888,549,4945  
6753 Engle Road  
Middleburg Heights, OH 44130  
[richmarweb.com](http://richmarweb.com)